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
GB – INSTRUCTIONS FOR USE

Description

The Coronet Donor Trephine Punch is a disposable punch and trephine system, for the trephination of the donor cornea button.

Intended Use

The Donor Trephine Punch components allow the surgeon to use a standard trephine to prepare and cut the donor button prior to transplantation. The graduated bowl of the base holds the donor cornea in an anatomically correct position while the trephine guide holds the trephine completely vertical during the procedure so creating a very clean vertical cut through the tissue.

CAUTION

- The device is for SINGLE USE ONLY. Do NOT re-sterilise or re-use.
- This device is intended for use by trained medical persons possessing the requisite skill and experience to use the device in accordance with the prevailing standards of medical practice and in conjunction with the instructions for this device.
- The punch is supplied sterile and ready for use. Do not use if the pouch has been opened or damaged.
- The Coronet trephine is much sharper than a traditional trephine. Caution should be exercised at all times when handling and using the trephine.
- The Coronet trephine is supplied with the blade in perfect condition. Should the device be accidentally dropped onto a hard surface, check the blade integrity under the microscope before use and discard the blade if damaged.
- The product comes into contact with bodily fluids, which can be contaminated. Care should be taken in the handling and disposal of the device after use to prevent contamination.
- CAUTION: US Federal law restricts this device to sale by or on the order of a physician.
- This product contains nickel which may cause allergic reaction in patients with sensitivity to nickel. Risk assess the use of the devices in relation to the medical benefit of the procedure and take necessary precautions with patients with known sensitivity to nickel.
- For safe disposal, the device should be considered as ‘sharp’ and must be immediately placed in a secure puncture-resistant bin (conforming to BS7320) which is suitable for incineration. All the components should be disposed of in accordance with hospital approved procedures.

Incident Reporting

- Any serious incident that has occurred in relation to this device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established

STERILISATION

- The device is a SINGLE USE ONLY device supplied sterile and ready for use. Sterilisation is by Ethylene Oxide (EO).
- DO NOT RE-STERILISE. Network Medical Products Ltd. assumes no liability for devices that have been re-sterilised by healthcare facilities.

SURGICAL TECHNIQUE

Separate the two parts of the punch - the base and the trephine guide.

Remove the cornea from the storage medium by grasping the rim of the surrounding sclera with toothed forceps and dry the anterior surface with a lint free PVA Sponge spear.

Position the cornea in the block well with the endothelial surface up. Centre the cornea using the guide holes.

Replace the trephine guide onto the base, being careful to maintain the centering of the cornea.

Remove the trephine from the packaging and position in the guide, taking care not to damage the ultra sharp cutting edge.

PRODUCT NAME DONOR TREPHINE PUNCH
PRODUCT INFORMATION AND INSTRUCTIONS FOR USE

Confirm the position of the corneal tissue and if necessary use a micro hook instrument via the windows of the trephine guide, to manipulate and position the corneal tissue optimally. Lower the trephine blade gently onto the endothelial surface, and then press with a gentle but firm movement on the trephine to achieve the vertical cut through the corneal tissue. It is not necessary to rotate the trephine but if in doubt the trephine and trephine guide can be rotated together to ensure a thorough cut.

Please note. The exceptionally super sharp trephine blade will cut through the corneal tissue without any sound or resistance. Once resistance is felt, the cut has been completed.

Once the cut has been achieved, remove the trephine and the trephine guide together in one piece. The corneal button, peripheral cornea and scleral ring will remain in the well of the base.

In the unlikely event that the corneal button remains in the lumen of the blade, a balanced salt solution or a viscoelastic substance should be gently dripped into the centre of the blade to dislodge the button.

HAZARDS ASSOCIATED WITH THE RE-USE OF SINGLE USE ONLY DEVICES:

1. Single use devices have not been validated for re-use.
If you re-use a device you may be held **Legally Liable** for the safe performance.
2. Cross-contamination and infection risks to patients . Including transmission of:
 - CJD & Variant CJD.
 - Prion Diseases.
 - Bacterial Endotoxins.
 - Hepatitis B & Hepatitis C.
 - Risks posed by HIV and AIDS
3. Device failure through material fatigue or degradation caused by initial use and design:
 - Plastics: Can be weakened, warped or become brittle.
 - Metals: Can be damaged or subjected to rusting. Blade edges will be damaged or blunt. Needles may be bent from the original shape.
 - Other materials: May degrade, becoming unacceptable when compared to original manufacturing criteria.
4. Patient injury from device failure and/or chemical burns from residue of decontamination agents absorbed into the materials.

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