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GB - INSTRUCTIONS FOR USE

Description

NETWORK External Nasal Splints are malleable padded aluminium splints which are cosmetically pleasing, easy to apply with hypoallergenic self-adhesive backing and require no heat to soften and mould to the shape of the patients' nose. Each splint is supplied with a dorsum pad to apply over the skin closure tape. Available in three sizes to accommodate for Caucasian, Asian and paediatric nasal anatomy.

Intended Use

An External Nasal Splint is for use after rhinoplasty or septoplasty procedures.

CAUTIONS

- This device is supplied NON STERILE and ready to use.
- This device is for SINGLE USE ONLY. Do NOT re-sterilise or re-use.
- Do not use if the packaging has been opened or damaged.
- This device is intended for use by trained medical persons possessing the requisite skill and experience to use the device in accordance with the prevailing standards of medical practice and in conjunction with the instructions for this device.
- The product comes into contact with bodily fluids, which can be contaminated. Care should be taken in the handling and disposal of the device after use to prevent contamination.
- CAUTION: US Federal law restricts this device to sale by or on the order of a physician.

Incident Reporting

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established

APPLICATION INSTRUCTIONS

1. IMPORTANT! Cleanse and dry skin on nose and application area.
2. Wipe cleansed area with Isopropyl Alcohol Wipe and dry thoroughly.
3. Express edema from application area.
4. Apply skin preparation to application area and DRY THOROUGHLY (Extremely IMPORTANT!)
5. Apply 1/2 inch wide skin closure tape starting at the bottom of the nose being careful to overlap each strip to be assured that all the skin in the application area is covered. (See Figure #1).
6. Remove release paper from Dorsum Pad and apply over skin closure tape from Radix to Supratip area if desired. Dorsum Pad adhesive should not contact skin. (See Figure #2)
7. Shape the Soft-Form AL splint to fit the application area.
8. Remove the release paper from the hypoallergenic adhesive on the back of the splint wings and apply foam side to nose over the Dorsum Pad and skin closure tape. Mould the Soft-Form Aluminium Splint to achieve proper fit and pressure. (See Figure #3)
9. Splint should be removed in seven to eight days.

Removal

1. Remove by cutting wings off with scissors, and gently lifting.
2. Retaping for three or four days is optional by the physician

Sterilisation

This device is supplied Non sterile and is designated for SINGLE USE ONLY.

HAZARDS ASSOCIATED WITH THE RE-USE OF SINGLE USE ONLY DEVICES

1. Single use devices have not been validated for re-use. If you re-use a device you may be held Legally Liable for the safe performance.
2. Cross-contamination and infection risks to patients. Including transmission of:
 - CJD & Variant CJD.
 - Prion Diseases.
 - Bacterial Endotoxins.
 - Hepatitis B & Hepatitis C.
 - Risks posed by HIV and AIDS
3. Device failure through material fatigue or degradation caused by initial use and design:
 - **Other materials:** May degrade, becoming unacceptable when compared to original manufacturing criteria.
4. Patient injury from device failure and/or chemical burns from residue of decontamination agents absorbed into the materials.

PRODUCT NAME: EXTERNAL NASAL SPLINTS - MALLEABLE
PRODUCT INFORMATION AND INSTRUCTIONS FOR USE

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