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GB - INSTRUCTIONS FOR USE

Description

NETWORK Thermoplastic External Nasal Splints are splints which are cosmetically pleasing, easy to apply with hypoallergenic self-adhesive backing. They require gentle heat in water to soften, after which they can be moulded to the shape of the patients' nose. Each splint is supplied with a dorsum pad to apply over the skin closure tape. Available in three sizes to accommodate for Caucasian, Asian and paediatric nasal anatomy.

Intended Use

An External Nasal Splint is for use after rhinoplasty or septoplasty procedures.

CAUTIONS

- This device is supplied NON STERILE and ready to use.
- This device is for SINGLE USE ONLY. Do NOT re-sterilise or re-use.
- Do not use if the packaging has been opened or damaged.
- This device is intended for use by trained medical persons possessing the requisite skill and experience to use the device in accordance with the prevailing standards of medical practice and in conjunction with the instructions for this device.
- The product comes into contact with bodily fluids, which can be contaminated. Care should be taken in the handling and disposal of the device after use to prevent contamination.
- CAUTION: Caution: US Federal law restricts this device to sale by or on the order of a physician.

Incident Reporting

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established

APPLICATION INSTRUCTIONS

Step 1: IMPORTANT: Cleanse and dry skin on nose and application area.

Step 2: Wipe cleansed area with Isopropyl Alcohol Wipe and dry thoroughly.

Step 3: Express edema from application area.

Step 4: Apply skin preparation to application area and DRY THOROUGHLY – THIS IS VERY IMPORTANT.

Step 5: Apply 1/2 inch wide skin closure tape starting at the bottom of the nose being careful to overlap each strip to be assured that all the skin in the application area is covered. (See Figure #1)

Step 6: Remove release paper from Dorsum Pad and apply over skin closure tape from Radix to Supratip area. Dorsum Pad adhesive should not contact skin. (See Figure #2)

Step 7: Size the splint to the application area. Splint may be trimmed with scissors if necessary. A properly sized splint should extend laterally to the lateral osteotomy sites, (but not beyond), from the Radix to the Supratip dome junction.

Step 8: Soften nasal splint in 65-75°C water (approximately 30 seconds).

Step 9: Remove the release paper from the back of the splint and apply adhesive side to nose over the Dorsum Pad and skin closure tape, hold in place until the splint cools and sets to form comfortably around the nose. Mould the Splint to achieve the proper fit and pressure. IMPORTANT NOTE: Do not let splint adhesive contact skin. (See Figure #3)

Step 10: Splint should be removed in seven to eight days.

Removal instructions

Step 1: Remove by gently spreading blunt bayonet type forceps between the skin closure tape and the patient's skin. Do not use excess force so that the skin tents up.

Step 2: Re-taping for three or four days is optional by the physician

Sterilisation

This device is supplied Non sterile and is designated for SINGLE USE ONLY.

HAZARDS ASSOCIATED WITH THE RE-USE OF SINGLE USE ONLY DEVICES

1. Single use devices have not been validated for re-use. If you re-use a device you may be held Legally Liable for the safe performance.
2. Cross-contamination and infection risks to patients. Including transmission of:
 - CJD & Variant CJD.
 - Prion Diseases.
 - Bacterial Endotoxins.
 - Hepatitis B & Hepatitis C.
 - Risks posed by HIV and AIDS
3. Device failure through material fatigue or degradation caused by initial use and design:

PRODUCT NAME: EXTERNAL NASAL SPLINTS - THERMOPLASTIC
PRODUCT INFORMATION AND INSTRUCTIONS FOR USE USEUSE

- **Other materials:** May degrade, becoming unacceptable when compared to original manufacturing criteria.
4. Patient injury from device failure and/or chemical burns from residue of decontamination agents absorbed into the materials.

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