

NETCELL® PVA Headbands



Product Overview

The NETCELL® PVA Headband has been specifically developed to provide an immediate cooling effect for the surgeon.

- Provide an immediate cooling effect via evaporation once hydrated with water
- Reduces contamination to the operating site by controlling and limiting perspiration
- Manufactured from ultra-smooth, biocompatible PVA sponge
- Soft, flexible and fibre free
- Supplied non-sterile, single use only

PVA Instrument Wipe Size Options

Format	Pack Size	Product Code
Standard	10	35-950

Material Specification

Product Component

Specification

Headband	100% Polyvinyl Alcohol Sponge (PVA), elastic
Pouch Packaging	100% virgin, high clarity polyethylene film
Outer Box/Carton Packaging	500 Micron Printed White Boxboard

Intended Use

The PVA headband is a single use disposable device, designed to provide an immediate cooling effect to the wearer once hydrated with water.

Instructions for Use

Instructions for use are not provided as this product can be used safely without.

Conformity to the European Directives

The range of PVA sponge products is supplied sterile and intended for short-term use in the nasal cavity. Therefore classification is Class I Sterile, (rule 5, Annex IX, 93/42/EEC Medical Devices Directive). Classification is therefore Class I STERILE