



Product Overview

The NETWORK ENT® House Adapter has been specifically designed to work with aspiration cannulas.

- Suction control
- Supplied sterile, single use only, one per pack, declared 5 year shelf life
- Manufactured in the UK

Product Size	Box Quantity	Product Code
House Adapter with Suction Control (Purple)	50	74-2001

House Adapter Specification Sheet

Material Specification

Product Component	Specification
Adapter	ABS
Pouch Packaging	Tyvek/Film
Outer Box/Carton Packaging	500 Micron Brown Boxboard

Intended Use

Handle with suction control hole used as a connector between tube and microsurgical tube to control suction manually and individually by the operator.

Sterilisation

Products are sterilised by Ethylene Oxide (EO) according to a validated cycle and the requirements of BS EN ISO 11135-1 (current).

Instructions for Use

The instructions for use and suggested surgical technique are supplied in the form of a multi-lingual leaflet with instructions given in diagram form where appropriate. The international symbols used on the packaging are explained in each language. A leaflet is supplied with each product.

Conformity to the European Directives

The Network range of Aspiration Suction Tubes and Fine Ends are identified as surgically invasive devices for transient use and are therefore classified as a Class IIa devices according to Annex IX, Rule 6, of the Medical Devices Directive 93/42/EEC as amended by 2007/47/EC.