




AN  
INNOVIA  
MEDICAL  
COMPANY




# Beaver<sup>®</sup> Otolaryngology Blades

Beaver<sup>®</sup> is part of the BVI Group of companies







# Beaver® Otolaryngology Products









## BEAVER® MINI-BLADES

	Description	Pack	Code
	<b>Beaver® 6400 Mini-Blade®</b> Grindless, rounded tip, sharp on one side, green	20 per box	<b>Beaver6400</b>
	<b>Beaver® 6600 Mini-Blade®</b> Grindless, lamellar, angled 60 degrees, smooth undersurface, sharp all around, bevel up, black	3 per box	<b>Beaver6600</b>
	<b>Beaver® 6610 Mini-Blade®</b> Grindless, lamellar, straight, smooth undersurface, sharp all around, bevel up, black	20 per box	<b>Beaver6610</b>
	<b>Beaver® 6700 Mini-Blade®</b> Grindless, smaller alternative to #15 general blade, orange	20 per box	<b>Beaver6700</b>
	<b>Beaver® 6900 Mini-Blade®</b> Grindless, straight, sharp all around, 180° cutting surface, double bevel blade, blue	20 per box	<b>Beaver6900</b>
	<b>Beaver® 6910 Mini-Blade®</b> Grindless, angled 60 degrees, sharp all around, 180° cutting surface, double bevel blade, blue	6 per box	<b>Beaver6910</b>

### Directions for Assembly

-  Hold the Beaver® Mini-Blade® tray between thumb and pointer finger.
-  Thread the Beaver® handle in a clockwise motion, onto the threaded tip.
-  Continue to thread handle until the thread is fully engaged and the Mini-Blade® is released from the tray clips.
-  Lift to remove Beaver® Mini-Blade® from tray.



## ENT/HEAD & NECK PROCEDURE BLADES

	Description	Pack	Code
	<b>Knife</b> Straight, sharp all around, 180° cutting surface double bevel blade	3 per box	<b>374769</b>
	<b>Needle Blade</b> 3.0mm cutting edge	3 per box	<b>375910</b>
	<b>Needle Blade</b> 5.0 mm cutting edge	3 per box	<b>375920</b>
	<b>Needle Blade</b> 7.0 mm cutting edge	3 per box	<b>375930</b>
	<b>Needle Blade</b> 1.5mm cutting edge	3 per box	<b>375940</b>
	<b>Beaver Mini-Blade</b> Sharp on two sides, pointed tip	6 per box	<b>376300</b>
	<b>Beaver Mini-Blade</b> Rounded narrow tip, sharp on one side	24 per box	<b>376401</b>
	<b>Tonsillectomy Blade</b> Sickle sharp	24 per box	<b>376800</b>

# Beaver® Otolaryngology Products

## BEAVER® HANDLES




Beaver Handles will accommodate all Beaver blades, supplied non-sterile.

	Description	Pack	Code
	Round, Knurled, Stainless 10.0 cm L x .50 cm D	1 per box	<b>374310</b>
	Round, Knurled, Stainless 13.0 cm L x .50 cm D	1 per box	<b>371350</b>

## MYRINGOTOMY BLADES AND KNIVES

	Description	Pack	Code
	<b>Myringotomy Blade</b> Lance tip, downward cut, offset 45°	6 per box	<b>377100</b>
	<b>Myringotomy Blade</b> Lance tip, upward cut, offset 45°	6 per box	<b>377110</b>
	<b>Myringotomy Blade</b> Spear tip, offset 45°, narrow shaft	6 per box	<b>377120</b>
	<b>Myringotomy Blade</b> Lance tip, downward cut, narrow shaft, offset 45°	6 per box	<b>377121</b>
	<b>Myringotomy Knife</b> Spear tip, offset 45°, narrow shaft	3 per box	<b>377122</b>
	<b>Myringotomy Knife</b> Lance tip, downward cut, narrow shaft, offset 45°	3 per box	<b>377124</b>

## TYMPANOPLASTY/STAPEDECTOMY BLADES

	Description	Pack	Code
	<b>Tympanoplasty Blade</b> 60° angled, sharp all around 2.5mm W, bevel down	3 per box	<b>377200</b>
	<b>Tympanoplasty Blade</b> Straight, sharp all around 2.5mm W	3 per box	<b>377210</b>
	<b>Sickle Edge Blade</b> Straight 84mm L 5.0mm cutting surface	3 per box	<b>377300</b>

All Beaver blades are supplied sterile



**Network Medical Products Ltd.**  
Coronet House, Kearsley Road,  
Ripon, North Yorkshire, HG4 2SG, UK.

Tel: +44 (0) 1765 609555 Fax: +44 (0) 1765 608476

email: [info@networkmedical.co.uk](mailto:info@networkmedical.co.uk)

[www.networkmedical.co.uk](http://www.networkmedical.co.uk)

Distributed by:

